


# Software

## Powerful microplate Reader Control and MARS Data Analysis Software


BMG LABTECH's software package includes Reader Control and MARS Data Analysis interfaces. This multi-user software package is included with every reader.

The intuitive Control Software is fully compliant with FDA regulation 21 CFR Part 11 and allows users to define instrument parameters and test protocols.

The MARS Data Analysis Software allows the user to display data, signal plots, spectra, and standard curves in 2D or 3D graphs. Data is processed using powerful predefined templates or a broad range of data calculation features. The software is also capable of creating standard curves and respective values (i.e.  $EC_{50}$ ,  $IC_{50}$  and  $r^2$ ) based on the following curve fitting algorithms:




Well scanning




NADH spectrum

- Linear regression
- 4- and 5-parameter
- Exponential
- Point-to-point
- Segmental regression
- Michaelis-Menten  $K_m$
- Cubic spline
- 2<sup>nd</sup> and 3<sup>rd</sup> polynomial.




Dose-response curve

The MARS wizard creates a step-by-step calculation of a standard curve, and important parameters such as S/N, Delta F % and Z' are easily obtained. Fast analysis of enzyme kinetic data using standard fit equations completes the MARS software package.



3D enzyme kinetic chart



Current state window

HTRF is a registered trademark of Cisbio Bioassays.  
DLR is a trademark of Promega Corporation.  
LanthaScreen is a registered trademark of Invitrogen Corp.  
Transcreeper is a registered trademark of BellBrook Labs.  
AlphaScreen, AlphaLISA, AlphaPlex, LANCE, and DELFIA are registered trademarks of PerkinElmer, Inc.  
© 2016 All rights reserved. All logos and trademarks are the property of BMG LABTECH.

# Applications

BMG LABTECH continuously works with all the leading reagent companies to optimize instrument settings for their existing assays and their newest chemistries.

BMG LABTECH's comprehensive searchable applications center reflects more than 25 years of expertise and innovations in microplate reading technology. Over 4000 references exemplify the flexibility and versatility of our readers, as well as their use in the chemical and biological sciences.

**Transcreeper®**  
FRET assays  
Binding studies  
DNA/RNA quantifications  
AlphaPlex™ technology  
HTRF®  
DLR™  
AlphaTechnology  
Protease activity  
Kinase activity  
SNP Genotyping  
LanthaScreen®  
HTS  
Apoptosis  
ORAC  
Enzyme activity  
Immunoprecipitation  
Dual luciferase assays  
Protein quantifications  
BRET assays  
Enzyme kinetics  
ROS detection  
Cell Viability  
PCR product quantifications  
NADH / NADPH assays  
DELFIA®  
Solubility tests  
LANCE®  
ATP and ADP detection  
ELISA  
Gene expression  
Reporter gene assays



**BMG LABTECH**  
*The Microplate Reader Company*

**Headquarters Germany**  
BMG LABTECH GmbH  
Allmendgrün 8  
77799 Ortenberg  
Tel. +49 781 96968 -0  
sales@bmglabtech.com

**Australia**  
BMG LABTECH Pty. Ltd.  
2/24 Carbine Way  
Mornington, Victoria 3931  
Tel. +61 3 5973 4744  
australia@bmglabtech.com

**France**  
BMG LABTECH SARL  
7, Rue Roland Martin  
94500 Champigny s/Marne  
Tel. +33 1 48 86 20 20  
france@bmglabtech.com

**Japan**  
BMG LABTECH JAPAN Ltd.  
2F TS-1 Building  
1-6-2, Shimo-cho  
Omiya-ku  
330-0844 Saitama City  
Tel. +81 48 647 7217  
japan@bmglabtech.com

**UK**  
BMG LABTECH Ltd.  
5 Alton House Office Park  
Gatehouse Way  
Aylesbury  
HP19 8XU  
Tel. +44 1296 336650  
uksales@bmglabtech.com

**USA**  
BMG LABTECH Inc.  
13000 Weston Parkway  
Suite 109  
Cary, NC 27513  
Tel. +1 877 264 5227  
usa@bmglabtech.com

www.bmglabtech.com

Made in Germany

# Instrumentation Guide




**BMG LABTECH**  
*The Microplate Reader Company*



# PHERAstar® FSX

## The New Gold Standard for HTS

Fluorescence Intensity  
 Fluorescence Polarization  
 Time-Resolved Fluorescence  
 TR-FRET/HTRF®  
 UV/Vis absorbance spectra  
 Luminescence  
 AlphaTechnology




3456 wells

### German engineering at its best.

The PHERAstar® FSX is the new Gold Standard microplate reader for High-Throughput Screening, specifically designed for highest sensitivity and speed. Its new and unique features make it superior to any other microplate reader currently on the market.

- Most sensitive reader in Fluorescence Intensity and Polarization
- Fastest read times with Simultaneous Dual Emission detection (incl. AlphaTechnology)
- 9 decades luminescence dynamic range
- Full absorbance spectra from 220-1000 nm in less than 1 s/well
- Top and bottom reading with focus adjustment (0.1 mm z-height)
- All microplate formats up to 3456-well
- New generation TRF laser for highest performance
- AlphaScreen®/ AlphaLISA® laser
- Up to two onboard reagent injectors
- Three integrated barcode readers



(R)-PFI 2 hydrochloride inhibitor titration with SET7/9 enzyme

### Never worry about which filter or dichroic mirror is installed!

Assay-specific Optic Modules are configured with all the necessary optical components including excitation and emission filters, dichroic mirrors, beam splitters, and polarization filters. The PHERAstar FSX can accommodate up to six Optic Modules. All Optic Modules are easily exchangeable, barcoded, and are automatically selected by the PHERAstar FSX for the appropriate assay.



# CLARIOstar®

## The LVF Monochromator™ reader

Fluorescence Intensity  
 Fluorescence Polarization  
 Time-Resolved Fluorescence  
 TR-FRET/HTRF®  
 UV/Vis absorbance spectra  
 Luminescence  
 AlphaTechnology




1536 wells

### The most sensitive monochromator-based reader

- Linear Variable Filters LVF Monochromators™
- Continuously adjustable bandwidths (8-100 nm)
- Increased sensitivity over conventional monochromators
- Fluorescence and luminescence spectral scanning
- Use monochromators, filters, or a combination of both
- Integrated fluorophore library for easy wavelength selection

### CLARIOstar® additional features

- Full absorbance spectra from 220-1000 nm in less than 1 s/well
- All microplate formats up to 1536-well
- Top and bottom reading with focus adjustment (0.1 mm z-height)
- Shaking and incubation up to 45°C or 65°C
- Up to two onboard reagent injectors
- AlphaScreen®/AlphaLISA® laser



Simplified representation of the CLARIOstar's LVF Monochromator.

### LVF Monochromator™ and filter selector technology

Highest light transmission and no stray light give the CLARIOstar's LVF Monochromator system filter-like performance and sensitivity.


The LVF Monochromator slides can also hold individual filters, polarizers, and dichroic mirrors for specific applications. This unique in-line optical path allows fixed filters to be combined with the LVF Monochromators in one measurement, providing users with the highest level of versatility, flexibility, and sensitivity.



# FLUOstar® Omega

## Filter based multi-mode reader

Fluorescence Intensity  
 Time-Resolved Fluorescence  
 TR-FRET/HTRF®  
 UV/Vis absorbance spectra  
 Luminescence  
 AlphaTechnology



384 wells

### Upgradeable multi-mode microplate reader

- Spectrometer-based absorbance with 220-1000 nm spectral scans in <1 sec/well, or filter-based absorbance with range 240-740 nm
- Top and bottom reading
- All microplate formats up to 384-well
- Plate shaking, incubation up to 45°C or 65°C
- Up to two onboard reagent injectors



# POLARstar® Omega

- Includes all FLUOstar Omega features
- Speed and performance increase with Simultaneous Dual Emission detection in FRET, BRET, and FP

# LUMIstar® Omega

- Dedicated luminometer for flash and glow assays
- Fully upgradeable to FLUOstar or POLARstar Omega
- Simultaneous Dual Emission for BRET



# NEPHELOstar® Plus

Unique laser-based nephelometer

### Light-scattering and turbidity measurements

- Uses light-scattering for detection of insoluble particles
- All microplate formats up to 384-well
- Up to two onboard reagent injectors
- Shaking and incubation up to 45°C
- Stacker and robot compatible



384 wells

# SPECTROstar® Nano

## Ultra-fast UV/Vis spectrometer

UV/Vis absorbance spectra  
 Microplate-based absorbance  
 Cuvette-based absorbance



1536 wells

### Ultra-fast DNA/RNA, protein, and ELISA measurements

- Full absorbance spectrum 220-1000 nm in <1 s/well
- All microplate formats up to 1536-well
- Well scanning, kinetic and endpoint measurements
- Integrated cuvette port, plate shaking and incubation up to 45°C
- Gas vent for atmospheric sensitive samples

## Accessories


### Atmospheric Control Unit (ACU)



#### Optimal environment for any live cell-based assay:

- Active regulation of O<sub>2</sub> and CO<sub>2</sub>
- Available for CLARIOstar

### LVIS Plate



#### Microplate for low-volume measurements

- Sixteen microdrop well sites for 2 µL samples
- Horizontal standard cuvette position
- NIST-traceable optical filters for precision and accuracy tests
- Compatible with all spectrometer-based readers

### Stacker



#### Mid-throughput microplate handling

- Rapid loading, unloading, restacking, for up to 50 microplates
- Continuous load feature and barcode reader
- Accommodates all microplate formats
- Compatible with CLARIOstar, PHERAstar FSX, Omega series, and NEPHELOstar Plus