

MACHEREY-NAGEL

Thin layer chromatography



Chromatography



State-of-the-art TLC products

- Efficient
- Reproducible
- Robust

Buy 5 – get 1 pack
815381.1 for free *



MACHEREY-NAGEL

www.mn-net.com



TLC promotion

What is TLC?

It's a quick, easy and cost-efficient method for separating compound mixtures utilizing planar plates or sheets made of glass, aluminum or polyester coated with silica or other separating media.

What is it used for?

A common method in pharmacy, food industries, natural compound analyses or organic syntheses for control of identity, yield or composition of substances and products.

Why should customers buy MACHEREY-NAGEL TLC products?

- Comprehensive range of plate sizes, surface chemistries and backings
- Pre-coated plates ready for immediate use
- Homogeneous, smooth and well adhering layers
- Consistent high quality from batch-to-batch and from plate-to-plate

Ordering information

Designation	REF	Plate size (cm)	Thickness of layer (mm)	Pack of
ALUGRAM® Xtra aluminum sheets				
ALUGRAM® Xtra SIL G/UV ₂₅₄	818329	2.5 x 7.5	0.2	200
ALUGRAM® Xtra SIL G/UV ₂₅₄	818331	4 x 8	0.2	50
ALUGRAM® Xtra SIL G/UV ₂₅₄	818332	5 x 20	0.2	50
ALUGRAM® Xtra SIL G/UV ₂₅₄	818333	20 x 20	0.2	25
Glass plates				
TLC plates SIL HD UV ₂₅₄	809220	10 x 10	0.25	25
TLC plates SIL HD UV ₂₅₄	809222	10 x 20	0.25	50
TLC plates SIL HD UV ₂₅₄	809223	20 x 20	0.25	25
TLC plates SIL HD UV ₂₅₄	809227	5 x 10	0.25	50
POLYGRAM® polyester sheets				
POLYGRAM® SIL G/UV ₂₅₄	805901	2.5 x 7.5	0.2	200
POLYGRAM® SIL G/UV ₂₅₄	805021	4 x 8	0.2	50
POLYGRAM® SIL G/UV ₂₅₄	805022	5 x 20	0.2	50
POLYGRAM® SIL G/UV ₂₅₄	805023	20 x 20	0.2	25

Promotion 2023

Buy 5 –
get 1 pack
815381.1
for free *

Details REF 815381.1

Standard silica gel
for liquid and flash
chromatography

(Silica gel 60 M, particle size:
0,04–0,063 mm, pack of 1000 g
in plastic container)



Ask for
better
price



Tech-Line SA · Chemicals-Quality-Control
4, Eleftherias Street · 14564 Kifisia · Greece
Tel.: +30 210 5281990-5 · Fax: +30-210-5222965



www.mn-net.com

MACHEREY-NAGEL



Management
System
EN ISO 13485:2016
ISO 9001:2015

www.tuv.com
ID 0000056401



MACHEREY-NAGEL GmbH & Co. KG · Valencienner Str. 11 · 52355 Düren · Germany

DE +49 24 21 969-0 info@mn-net.com
CH +41 62 388 55 00 sales-ch@mn-net.com

FR +33 388 68 22 68 sales-fr@mn-net.com
US +1 888 321 62 24 sales-us@mn-net.com